

GF S.p.A. News Letter




Dear sirs,

GF S.p.A. is pleased to present one of our newest inspection machines for prefilled syringes. This represents a rotary inspection unit with continuous automatic movement, for single-dose syringes. By means of special grippers, the syringes are extracted from the infeed star and inserted in the carousel where they will undergo inspection checks. After they have been inspected, the syringes are placed on the outfeed star and are then ready to be directed into the various outfeed channels.



Type of checks carried out

1. Particle inspection
2. Cosmetic inspection
3. Fill level check
4. Plunger inspection
5. Tip cap inspection
6. Inspection of syringe bottom integrity and fingergrip



Please do not hesitate to contact us for any further information.

Best regards



Pharmaceutical Machinery Division

G.F. S.p.A. - Via T. Edison, 3 - 43045 Rubbiano fraz. Solignano (Parma) - Italy

Tel. +39 0525 300311 - Fax +39 0525 300350 - Internet: <http://www.gf-industries.it> - E-mail: gfe@gfe.it